



Z2® Accessories

The following Accessories are available at your retail medical supplier.



PowerShell® With Extended Life Battery (005767) - The PowerShell integrates the Z2 with a rechargeable lithium ion battery, so you can take your CPAP anywhere, even where there is no power.




Extended Life Battery (006910) - The Extended Life Battery is a rechargeable, lithium ion battery that integrates with the PowerShell to deliver battery power to the Z2.




Z2 Bag (009001) - Protect your Z2 with this padded bag with a drawstring closure. Lays flat for easy storage when not in use.



HME (Heat Moisture Exchange) (005998) - A small, portable, disposable humidification cartridge.



4 ft Slim Style Tube (005797) - 4 foot lightweight, thin pressure delivery tube for use with the Z2 device.



6 ft Slim Style Tube (006908) - 6 foot lightweight, thin pressure delivery tube for use with the Z2 device.




Q-lite® In-line CPAP Muffler (007961) – Reduces acoustic noise transmitted to the patient mask through the CPAP tube. Includes an 8” extension tube.



Qtube® Replacement Kits (006005) – Contains the replacement foam, O ring, and plastic wrap for the original Qtube Muffler.



Filters (2 pack) (005788) – Replacement polyester air inlet filters for the Z2.



Tube Adapter (005779) – Required for Z2 operation. Adapts a standard 22mm CPAP hose to the outlet of the Z2.



USB Cable (005757) – USB to micro-USB data cable enables data transfer between a PC and Z2.



Plastic End Cap (005775) – Replacement end cap fits the Z2.



Power Supply (005755) – 100-240V power supply for the Z2.

Please note, these Accessories can be used with either the Z2 or the Z1® CPAP from Breas Medical. For more information visit www.breas.us.

Z2®, Q-lite®, PowerShell® and Qtube® are registered trademarks of Breas Medical, Inc. Breas and the Breas logo are trademarks of Breas Medical, AB. © 2024 Breas Medical - All rights reserved. BOS-SAL-MAR-000333 Rev4

Breas Medical Inc · 16 Esquire Rd · N. Billerica · MA 01862 · USA · Phone 855-436-8724 · www.breas.us