




EveryWare – Manage your respiratory patients remotely with confidence

EveryWare is a securely-hosted cloud-based application offering clinicians & care givers EveryWare-access to patient therapy data.



By enabling the provision of care outside of the hospital, the aim of EveryWare is to help make a more productive use of healthcare resources and reduce expensive secondary care admissions.

SECURELY HOSTED

EveryWare is a securely-hosted cloud-based application offering clinicians & care givers EveryWare-access to patient therapy data.

4G CELLULAR CONNECTION

Transmitted using a 4G cellular connection, Breas patients can be remotely managed with confidence by demonstrating their device usage whilst meeting the latest cyber-security compliance needs.

MANAGEMENT BY EXCEPTION

Built on a management-by-exception approach, EveryWare's intuitive user interface and presentation of therapy data aims to generate intelligence-led insights, which in turn can help optimise workflows and ultimately deliver better patient outcomes.

At Breas, we know that ageing populations and modern lifestyles can lead to chronic conditions that stretch healthcare access, resources and budgets to crisis point.

Chronic respiratory conditions can present patients with complex needs who require a multi-disciplinary team to manage them effectively.

Remote technology has the opportunity to not only enhance patient's quality of life and make their treatment more effective, but also help plan interactions with patients more effectively, cutting out needless and inconvenient travelling.



For tutorial videos and other educational material on ventilation and airway clearance, please scan or click the QR code to visit Education by Breas.



Introducing change settings, the latest update to EveryWare by Breas.

Changing Settings for Breas' non-life support ventilators exemplifies EveryWare's core value of improving Outcomes through Optimisation.

Recognising that each patient will respond differently to treatment, caregivers may need adjust prescription settings to help deliver maximum comfort and optimise treatment outcome for their patients. To date, these adjustments often necessitated multiple face-to-face appointments with the ventilator in front of them.

This process can now be accelerated thanks to the ability to Change Settings remotely, direct from within EveryWare. The time previously spent can now be used for managing more complex patients, or to make a more productive use of healthcare resources.

EveryWare now allows caregivers to remotely change settings ranging from IPAP/EPAP to Inspiratory Time (Ti) and High Flow Nasal Therapy (HFNT) flow* and alarms ranging from Low Pressure through to High EtCO₂.



Vivo 1



Vivo 2



Vivo 3

*Vivo 3 only

Change Settings Requirements:

- EveryWare 2.0.x with an authorised User Role.
- Compatible Breas Non-Life Support ventilator and firmware.
- Connected to iLink with cellular connection.