



## The PowerShell® with Extended Life Battery

The PowerShell is the first fully integrated power solution for CPAP users that need a totally portable, integrated power source. It is small enough to fit in the palm of your hand but has enough power for a full night of cord-free CPAP. The Extended Life Battery should provide 8.9 hours of CPAP therapy at the highest settings (20 cmH2O, Z-Breathe® 3). The PowerShell is compatible with both the Z1® and Z2® CPAPs.



- Soft, padded 2mm neoprene skin over a solid plastic skeleton ensures a secure fit
- Surrogate power plug for built-in charging using the standard Z1/Z2 power adapter
- Insert the Extended Life Battery module and the Z2 or Z1 to create an integrated form-factor

### Breas Extended Life Battery Technical Specifications

Specifications	
Weight	15.5oz / 440g
Power	14.4 V
Capacity	99 Watt Hours (Wh)
Operating Temp	41°-95°F/5°-35°C



#### Intended use:

Your Z1 or Z2 Auto System is a single patient, reusable device that provides continuous positive airway pressure (CPAP) to support treatment of adults weighing over 66 lbs (30 kg) with obstructive sleep apnea (OSA).

Z2®, Z1®, Z-Breathe® and PowerShell® are trademarks of Breas Medical Inc. © 2024 Breas Medical – All rights reserved. Breas Medical AB reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. Breas and the Breas logo are trademarks of Breas Medical AB. COR-SAL-MAR-008066-Rev0