



Product Update: New Trolley Now Available

New! GCX Vivo 45LS Trolley (PN 010472) replaces Breas Ventilator Trolley (007384)

Introducing our new trolley, designed for enhanced durability, smooth, stable mobility, and easy assembly.

Breas' new trolley: GCX Trolley (PN 010472)

- Durable, larger, high-quality 360° caster wheels enable smooth and effortless mobility
- Simplified assembly
- Validated for use with the Vivo 45LS
- Load Capacity: Maximum externally added load: 12.1 kg (26.6 lb) in standard configuration
- Existing compatible accessories: Vivo 45LS mounting bracket, IV-pole, circuit support arm
- Product of the USA

Discontinued: Breas Trolley (PN 007384)

- Discontinued as of November 15, 2025

Replacement Parts and Accessory Compatibility

PN	Description	Status
010472	Trolley GCX for Vivo 45 Series	New
010590	Screw kit for Trolley GCX (010472)	New
006761	Vivo 45LS mounting bracket	Compatible
007859	Trolley IV-pole	Compatible
007917	Trolley Circuit Support Arm	Compatible
007858	Universal Rail Clamp (Trolley)	Compatible
005128	Oxygen Cylinder Holder (Trolley)	Compatible



New GCX Trolley available from Breas.

BREAS.US

Not compatible with Vivo 2, Vivo 3, Vivo 65, or Clearo.

Vivo, Vivo by Breas, Breas, and the Breas logo are trademarks of Breas Medical AB. © 2025 Breas Medical, Inc. All rights reserved. Breas Medical AB reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. MAR-008189 Rev 0

Breas Medical Inc · 16 Esquire Rd · N. Billerica · MA 01862 · USA · Phone 855-436-8724 · www.breas.us