



## Nitelog – Simplicity at your fingertips

Nitelog® wirelessly syncs Z2® Auto CPAP sleep data to an iOS or Android device to analyze or share. The new and improved interface is easy to use. It makes data sharing with healthcare providers easy.

Nitelog connects to your Z2 Auto CPAP using the built-in wireless connectivity. Your sleep data can be saved to the cloud to back up your data and restore it to another wireless device.



### ANALYZE

It's never been easier to oversee and interact with your sleep data. Check your progress from anywhere right in the palm of your hand. Set your sleep goals and get feedback on your sleep quality with colored AH index emojis.



### SHARE

Nitelog is a simple, easy way to share sleep data with providers. Just generate a compliance report and email it directly from your phone. Z1® devices may need to be updated to be compatible.



### IPHONE AND IPAD USERS

To download Nitelog on an iOS device, go to the App Store on that device (iPhone or iPad) and search for "Nitelog". Push "GET" to download.




### ANDROID USERS

To download Nitelog on an Android device, go to the Google Play store on that device (Android phone or tablet) and search for "Nitelog". Press the "INSTALL" button to download.

# Breas Nitelog Technical Specifications

Specifications	
Platforms	Apple / Android
Connecting Devices	Z1, Z2 Auto version
Operating Modes	Fixed pressure CPAP, Auto-adjusting APAP
Settings	Z-Breathe®, Auto Start, APAP or CPAP Pressure
User Settings	Ramp Time, Ramp Start Pressure, Sleep Start Time, Sleep Goals
Compliance Report Contents	Usage Date range, Required Hours per Day, Required Percentage of Compliant Days, Days in Range, Number of Days Used, Number of Compliant Days, Percentage of Compliant Days, Compliant Pass/Fail, Pressure Average, AHI, AHI Average, Mode, Minimum Pressure, Maximum Pressure, CPAP Pressure, Ramp Time, Ramp Start Pressure, Z-Breathe® Setting

## Screens



### Intended use:

Your Z1 or Z2 Auto System is a single patient, reusable device that provides continuous positive airway pressure (CPAP) to support treatment of adults weighing over 66 lbs (30 kg) with obstructive sleep apnea (OSA).

Z2®, Z1®, Z-Breathe® and Nitelog® are trademarks of Breas Medical AB. © 2024 Breas Medical – All rights reserved. Breas Medical reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. Breas and the Breas logo are trademarks of Breas Medical AB. BOS-MAR-MAR-008060 Rev0