

BREAS

UK Breathing Circuits
and Consumables



Breas UK - Breathing Circuits Feedback Form




Breas need to capture market feedback on all of its products. If you use a Breas breathing circuit, please take a moment to complete this short questionnaire to provide us with your experiences and thoughts.

Just scan or click the above QR code to provide your feedback.

Page	Contents
1.	<u>Invasive 22mm Circuits - Adult</u>
2.	<u>Invasive 22mm Circuits - Paediatrics</u>
3.	<u>Invasive 15mm Circuits - Paediatrics</u>
4.	<u>Non-Invasive 22mm Circuits</u>
5.	<u>Non-Invasive 15mm Circuits</u>
6.	<u>22mm Non-Invasive Circuit (acute use)</u>
7.	<u>Invasive 22mm Active Valve Circuits</u>
8.	<u>Miscellaneous Ventilation Circuits</u>
9.	<u>Ventilation Circuits - Additional Items</u>
10.	<u>Mouthpiece Ventilation Circuits/Kits</u>
11.	<u>15mm and 22mm CPAP tubing</u>
12.	<u>Nippy Clearway and Clearway 2 circuits</u>
13.	<u>Clearway interfaces and consumables</u>


Invasive 22mm Humidified/Heated Leakage Circuit for use with MR850 Humidifier with Adult Catheter Mount

Product Code: 008056 | Box Quantity of 10



Invasive 22mm Dry Leakage Circuit with Adult HME and Catheter Mount


Product Code: 008050 | Box Quantity of 10



[BACK TO CONTENTS PAGE](#)


Invasive 22mm Humidified/Heated Leakage Circuit for use with MR850 Humidifier with Paediatric Catheter Mount

Product Code: 008059 | Box Quantity of 10



Invasive 22mm Dry Leakage Circuit with Paediatric HME and Catheter Mount


Product Code: 008048 | Box Quantity of 10



[BACK TO CONTENTS PAGE](#)

Invasive 15mm Humidified/Heated Leakage Circuit for use with MR850 Humidifier


Product Code: 009696 | Box Quantity of 10




- | | | | | | |
|---|---|---|---|---|--------------------------------------|
| A | - | 1690001 Anti Bacterial Filter | E | - | Humidifier Temperature Monitor Point |
| B | - | 5006000 22mm Smoothbore Tubing (600mm length) | F | - | 5802001 Purple CO2 leak port |
| C | - | Humidifier Temperature Monitor Point and Heated Wire Connection | G | - | Swivel Catheter Mount |
| D | - | 5015008 15mm Smoothbore heated Tubing (1500mm length) | H | - | Auto Fill Water Chamber |

Invasive 15mm Dry Leakage Circuit

Product Code: 008082 | Box Quantity of 10




- | | | |
|---|---|--|
| A | - | 1690001 Intersurgical Bacterial Filter |
| B | - | 5016000 Intersurgical 15mm Smoothbore tubing |
| C | - | 1962000 Intersurgical 15 -22F connector |
| D | - | 5802001 Intersurgical Exhalation Port 22M -15F/22M |
| E | - | 1850000 Intersurgical HME 15M 15F/22M |
| F | - | 3535000 Intersurgical Catheter Mount |




Non-Invasive 22mm Dry Circuit for use with Vented Interface

Product Code: 008024 | Box Quantity of 10



Non-Invasive 22mm Humidified Circuit for use with Vented Interface


Product Code: 008111 | Box Quantity of 10



[BACK TO CONTENTS PAGE](#)


Non-Invasive 15mm Dry Circuit for use with Vented Interface

Product Code: 008028 | Box Quantity of 10



Non-Invasive 15mm Humidified Circuit for use with Vented Interface


Product Code: 008113 | Box Quantity of 10



[BACK TO CONTENTS PAGE](#)


Non-Invasive 22mm Dry Leakage Circuit for use with Non-Vented Interface

Product Code: 008026 | Box Quantity of 10



Non-Invasive 22mm Humidified/Heated Circuit for use with Vented Interface


Product Code: 008074 | Box Quantity of 10



[BACK TO CONTENTS PAGE](#)


Invasive 22mm Humidified/Heated IPPV Circuit for use with MR850 Humidifier

Product Code: 008101 | Box Quantity of 10



Invasive 22mm Dry IPPV Circuit

Product Code: 008034 | Box Quantity of 10



Pilot Line Connector for using Nippy 4+ with an active exhalation valve circuit

Product Code: 007654

Supplied with every Nippy 4+, or sold individually

[BACK TO CONTENTS PAGE](#)

Miscellaneous Ventilation Circuits



22mm SeptiFlow Patient Circuit with active exhalation valve

Product Code: 007616

Sold individually



Pilot Line Connector for using Nippy 4+ with an active exhalation valve circuit

Product Code: 007654

Sold individually



Heated Patient Circuit (Vivo 1-2-3)

Product Code: 006990

Sold individually



Heated Patient Circuit (Nippy 4 / 4+)

Product Code: 006193

Sold individually

[BACK TO CONTENTS PAGE](#)

Ventilation Circuits - Additional Items



Anti Bacterial Filter

Code: 008018

Sold individually



Swivel Leakage Port

Code: 004426


Sold individually



HC325 Manual Fill Water Chamber

Code: 008351

Box Quantity of 10



Extendable Paediatric Catheter Mount (with suction port)

Code: 008307


Sold individually



15mm Heated Moisture Exchange (HME)

Code: 008019

Sold individually



Extendable Swivel Catheter Mount (with suction port)

Code: 008306


Sold individually



22mm Heated Moisture Exchange (HME)

Code: 008021

Sold individually



Extendable Swivel Catheter Mount (with suction port)

Code: 008302


Sold individually



Leakage Port

Code: 008130

Sold individually



Extendable Swivel Catheter Mount (no suction port)

Code: 008303

Sold individually

[BACK TO CONTENTS PAGE](#)

Mouthpiece Ventilation Kits - UK

The Nippy 4 can deliver an MPV-Pressure mode, and the Nippy 4+, Vivo 55 and Vivo 65 ventilators are all capable of delivering MPV-Pressure and Volume modes for your patients. You can find instructions in the User Manuals, and tutorial videos explaining the use of MPV can be found by clicking on or scanning the QR code. Breas have all of the required accessories available for patients using Mouthpiece Ventilation.



SCAN ME




Mouth Piece Ventilation Circuit
fully assembled



008254

MPV Flexible Arm Assembly



006094 - Mouthpiece (15 pcs)

006093 - MPV Circuit, Disposable (2 pcs)

CPAP tubing



22mm Lightweight Smoothbore Breathing Circuit

Product Code: 008426

Not sold individually - only available in box quantity of 30

Discounts available for agreed annual volumes



15mm Lightweight Smoothbore Breathing Circuit

Product Code: 008427

Not sold individually - only available in box quantity of 30


Discounts available for agreed annual volumes

[BACK TO CONTENTS PAGE](#)

Nippy Clearway Circuits


MI-E circuit

Product Code: 008157 | Box Quantity of 10



IPPB circuit


Product Code: 008185 | Box Quantity of 10



Clearway 2 Circuits


MI-E circuit

Product Code: 008276 | Box Quantity of 10



IPPB/NIV circuit

Product Code: 008277 | Box Quantity of 10



[BACK TO CONTENTS PAGE](#)

Nippy Clearway and Clearway 2 Interfaces and consumables



Mask Size 6 (XL Adult)

Product Code: 008342

Sold Individually



Mask Size 5 (Large Adult)

Product Code: 008341

Sold Individually



Mask Size 4 (Adult)

Product Code: 008340

Sold Individually



Mask Size 3 (Small Adult)

Product Code: 008339

Sold Individually



Mask Size 2 (Paediatric)

Product Code: 008338

Sold Individually



Mask Size 1 (Infant)

Product Code: 008337

Sold Individually



Mask Size 0 (Neonate)

Product Code: 008336

Sold Individually



Mouthpiece

Product Code: 008160

Box Quantity of 10



Mouthpiece and Nasal Clip

Product Code: 008161

Box Quantity of 10



Catheter Mount

Product Code: 008159

Box Quantity of 10



Anti Bacterial Filter

Product Code: 008169

Box Quantity of 10



Oxygen Coupler

Product Code: 008170

Box Quantity of 10



CW2 In-let Filter, White

Product Code: 009065

Pack of 5



CW2 In-let Filter, Grey

Product Code: 009066

Pack of 5

Breas UK - Breathing Circuits

Feedback Form



© 2023 Breas Medical Ltd – All rights reserved. Breas Medical Ltd reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. Breas and the Breas logo are trademarks of Breas Medical AB.

Breas Medical Ltd. Unit A2 The Bridge Business Centre, Timothy's Bridge Road, Stratford upon Avon, Warwickshire. CV37 9HW. United Kingdom. +44(0)1789 293 460 Company registered in England and Wales No. 07190611.

Document Number STR-MAR-MAR-008019-Rev_3