

# Vivo 2 - Adaptable Bi-Level

The Vivo 2 is a bi-level ventilator designed to prioritize patient comfort and ease of use. It offers personalized and comfortable non-invasive respiratory support for non-dependent patients with chronic breathing insufficiency. The intuitive interface ensures straightforward operation, making it user-friendly and accessible. In addition, a comprehensive set of modes and extensive monitoring capabilities enhances communication between the device and healthcare professionals, allowing for greater flexibility and ensuring optimal patient care.



### COMPACT

- One platform three models
- Lightweight and compact
- Integrated humidifier with heated wire circuit



#### **ADAPTABLE**

- Wide range of modes
- Adult and paediatric (>10 kg)
- Intuitive interface
- Extensive therapy data
- Suitable for use in hospital and home
- Closer monitoring with Breas accessory box



### CONNECTED

Making connected care more insightful

- EveryWare cecure cloud-based remote connectivity solution
- PC Software desktop-based with Live Viewer



## Breas Vivo 2 Technical Specifications

Settings/Performance	
Ventilation Modes	S, ST, T, PSV, PCV+A , CPAP
Settings Profiles	Three profiles
IPAP	4-30 cmH <sub>2</sub> O
EPAP	2-20 cmH <sub>2</sub> O
Breath Rate/Back-up Rate	0-40 bpm
Inspiratory Time	0.3-3 s
Rise Time	1-9, Auto is similar to setting 3
Inspiratory Trigger	1-9, Auto is similar to setting 3
Expiratory Trigger	1-9, Auto is similar to setting 5
Minimum Inspiratory Time	Off, 0.3-3 s
Maximum Inspiratory Time	Off, 0.3-3 s
Target Volume	Off, 100-1500 ml
Ramp (Up/Down)	Off, 10-60 min
Auto-EPAP	On/Off, Pressure Limit: Max EPAP +2 cmH <sub>2</sub> O-30 cmH <sub>2</sub> O, PS 2-28 cmH <sub>2</sub> O, Min/Max EPAP: 2-20 cmH <sub>2</sub> O, EPAP Step: 0.5-2 cmH <sub>2</sub> O, Relax Time: 2-12 min/Off
Humidifier Level	1-5
Heated Circuit	Off, 1-5
Patient Circuit	15 mm, 22 mm
O <sub>2</sub> Bleed in (Device)	0-30 l/min
Monitoring	
Displayed Data	Ppeak, EPAP, Leakage, MVe, Vte, % in TgV, Total Rate I:E, Spont Rate, % Spont, Inspiratory Time, Rise Time, SpO <sub>2</sub> , Pulse, Rate
Power Supplies	
Mains Supply	100-240 V AC
External DC	12-24 V
Xpac Extended-life Backup Battery (optional)	36 hours (IPAP 20 cm/EPAP 8 cm, Rate 15 bpm)
Alarms	
Alarms	High Pressure, Low Pressure, High Breath Rate, Low Breath Rate, High MV, Low MV, Rebreathing, Apnea, Disconnection, High EPAP, Low EPAP, High SpO <sub>2</sub> , Low SpO <sub>2</sub> , High Pulse Rate, Low Pulse Rate, Power Failure, High Pressure Limitation
Dimensions	
$W \times H \times D$	166 x 185 x 200 mm
Weight	1.6 kg
Noise Level (at 10 cmH <sub>2</sub> O constant pressure)	Less than 28 dB(A)
Connectivity and Software	
Built-in Connectivity	USB
Compatible Software	PC Software, EveryWare Remote Connectivity Solution
Plug-and-Play Accessories	SpO <sub>2</sub> , Nurse Call, Remote Alarm, Accessories Box
Other	
Humidification	Integrable water chamber
Circuits	Cleanable & disposable leakage circuits with/without heated wires
Confirm Stop	On, Off
Compliance/Therapy Data (on Screen)	Compliance/Usage Data, AHI, AI, HI, IPAP, EPAP, Leakage, Vte, MVe, Total Rate, % Spont Breaths, Insp Time, I:E

#### Intended use

The Vivo 2 is intended to provide non-invasive ventilation for patients weighing over 10 kg (22 lbs) who require long-term support or mechanical ventilation for respiratory impairment, with or without obstructive sleep apnea. The Vivo 2 is intended for spontaneously breathing patients.

Breas and the Breas logo are trademarks of Breas Medical AB. Vivo by Breas is a trademark of Breas Medical AB. © 2024 Breas Medical AB - All rights reserved. Breas Medical AB reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. COR-SAL-MAR-000363-Rev7