



Which Z1 is right for you?

Z1 CPAP

Z1 Auto

Fixed Pressure Mode	●	●
Auto Adjusting Pressure Mode		●
Auto Stop Feature	●	
Auto Start/Stop Feature		●
Approved for in-flight use	●	●
Mobile Data Capable		●
Pressure Range 4-20 cmH <sub>2</sub> O	●	●
Ramp 0-45 minutes	●	●
Z-Breathe® Comfort Feature	●	●
Q-Lite optional	●	●
No special mask required	●	●
Auto Adjusting Voltage for world travel (100-240v)	●	●
Weight	292g	299g

Dimensions: 16.5 x 8.4 x 5.1 cm

**The PowerShell® - Integrated Battery Power to-go.**  
The optional PowerShell® accessory is an integrated cpap battery solution that delivers a full night's sleep on a single charge. No cords required.



### Modular Batteries

Need more nights of off-the-grid power? Buy an additional battery module.



#### Battery Specifications

Battery capacity	99.4 Wh
Operating temperature	5°-35°C (41°-95°F)

### Explore your sleep data with the Nitelog® App

With the Nitelog 2.0 app you have all of your sleep data available at your fingertips. The improved interface is easy to use and connects to your Z1 Auto using the built-in Bluetooth chip. Nitelog also enables sharing by email, just generate a compliance report and email it directly to your doctor.



©2023 Breas Medical - All rights reserved. Breas Medical reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. Breas and the Breas logo are trademarks of Breas Medical AB. COR-SAL-MAR-001112-Rev7

Breas Medical AB · Företagsvägen 1 · SE-435 33 Mölnlycke · Sweden  
Phone +46 31 86 88 00 · breas.com



# Z1 Auto

Ultra-small, travel CPAP



# Take your adventures off the grid

Whether you need the fixed pressure Z1 CPAP or the auto-adjusting Z1 Auto, you can take your CPAP anywhere you want to go.

The optional PowerShell® with Extended Life Battery provides battery power for a full night's sleep on a single charge. It's great for camping or any where access to power is limited.

## Travel anywhere in the world

Both Z1s are FAA Compliant for in-flight use. They come with auto-adjusting voltage power supplies that let you plug into any outlet in the world with just an added plug adapter.

