

## Z2 Auto


The Z2 Auto CPAP machine is ultra-small, light weight, full featured auto pressure CPAP machine that is perfect for travel. The therapy pressure adjusts on a breath-by-breath basis which results in a lower overall mean pressure than a fixed pressure CPAP. It is the ultimate CPAP for travelers and people with an active lifestyle.

### Z2 Auto

Fixed Pressure mode	●
Auto Adjusting Pressure mode	●
Auto Start/Stop Feature	●
Approved for in-flight use	●
PC data capable	●
Mobile data capable	●
Pressure Range 4-20 cm H <sub>2</sub> O	●
Ramp 0-45 minutes	●
Z-Breathe® Comfort Feature	●
Auto Adjusting Voltage for world travel (100-240v)	●
Qtube included	●
No special mask required	●

## Check your treatment information

The backlit LCD Screen shows your treatment settings on the device. You can easily adjust your comfort settings directly on your device.



1. Treatment pressure display
2. Ramp progress indicator (if ramp set)
3. Ramp time remaining (if ramp set)
4. Battery status (if present)
5. AC Power indicator (if connected)
6. General status display (date, mode, Auto Start/Stop)
7. Fault indicator (if fault occurs)

**If you travel for work or just for fun  
ask about the Z2 today.**

©2023 Breas Medical - Z2™ and Z2 Auto™ are trademarks of Breas Medical, Inc. Nitelog® and Z-Breathe® are registered trademarks of Breas Medical, Inc. N. Billerica, MA. Breas and the Breas logo are trademarks of Breas Medical, AB. BOS-SAL-MAR-000337-Rev2

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., Google Play is a trademark of Google Inc., Bluetooth is a registered trademark of Bluetooth SIG, Inc.

Breas Medical, Inc. · 16 Esquire Road · North Billerica · MA01862 · USA  
Phone 855 436 8724 · [www.breas.us](http://www.breas.us)



## Z2 Auto

Ultra-small, travel CPAP



**Quiet motor  
Auto Start/Stop  
No special mask required  
Battery option**



## Take your adventures off the grid

With the auto-adjusting Z2 Auto, you can take your CPAP anywhere you want to go. Just add the optional PowerShell® with Extended Life Battery to provide battery power for a full night's sleep on a single charge.

It's great for camping or any where access to power is limited.

Shown actual size  
6.48" x 3.30" x 2.02"



## Travel anywhere in the world

The Z2 is FAA compliant for in-flight use. It comes with an auto-adjusting voltage power supply that lets you plug into any outlet in the world with just an added plug adapter.

## Check your sleep data

Download your Z2 Auto® sleep data directly to your iOS or Android phone or tablet with the free Nitelog® mobile app.

Just download the app, connect to the Z2 via Bluetooth and you're ready to go.

